


Biohit Proline Mechanical

Biohit Proline mechanical pipettors deliver high quality with excellent value for money. The Proline range is manufactured from the latest materials ensuring a light weight and durable pipettor which will give many years of reliable service. The Proline pipettors cover the volume range 0.1 µl to 5000 µl with 8 different single-channel and 6 multichannel adjustable volume models. The Proline is also available in 12 fixed volume models.



Proline Mechanical Offers

- Lightweight ergonomic design
- Easy calibration and maintenance
- Full documentation, warranties and certificates
- Colour-coded ID stickers included with each pipettor



Proline Single- and Multichannel, Performance Specifications and Ordering Details

Cat. No.	Channels	Volume Range (µl)	Increment (µl)	Test Volume (µl)	Inacc. (%)	Impr. (%)	Safe-Cone filters 50 pcs	
							Standard	Plus
720005	1-channel	0.1-2.5 µl	0.05 µl	2.5 1.25 0.25	2.50 % 3.00 % 12.00 %	2.00 % 3.00 % 6.00 %	-	-
720000	1-channel	0.5-10 µl	0.1 µl	10 5 1	1.00 % 1.50 % 2.50 %	0.80 % 1.50 % 1.50 %	-	-
720080	1-channel	2-20 µl	0.5 µl	20 10 2	0.90 % 1.20 % 3.00 %	0.40 % 1.00 % 2.00 %	721008	721018
720020	1-channel	5-50 µl	0.5 µl	50 25 5	0.60 % 0.90 % 2.00 %	0.30 % 0.60 % 2.00 %	721008	721018
720050	1-channel	10-100 µl	1 µl	100 50 10	0.80 % 1.00 % 3.00 %	0.15 % 0.40 % 1.50 %	721007	721017
720070	1-channel	20-200 µl	1 µl	200 100 20	0.60 % 0.80 % 3.00 %	0.15 % 0.30 % 1.00 %	721007	721017
720060	1-channel	100-1000 µl	5 µl	1000 500 100	0.60 % 0.70 % 2.00 %	0.20 % 0.25 % 0.70 %	721006	721016
720110	1-channel	1000-5000 µl	50 µl	5000 2500 1000	0.50 % 0.60 % 0.70 %	0.15 % 0.30 % 0.30 %	721005	721015
720210	8-channel	0.5-10 µl	0.1 µl	10 5 1	1.50 % 2.50 % 4.00 %	1.50 % 2.50 % 4.00 %	-	-
720220	8-channel	5-50 µl	0.5 µl	50 25 5	1.00 % 1.50 % 3.00 %	0.50 % 1.00 % 2.00 %	721014	-
720240	8-channel	50-300 µl	5 µl	300 150 50	0.70 % 1.00 % 1.50 %	0.25 % 0.50 % 0.80 %	721014	-
720310	12-channel	0.5-10 µl	0.1 µl	10 5 1	1.50 % 2.50 % 4.00 %	1.50 % 2.50 % 4.00 %	-	-
720320	12-channel	5-50 µl	0.5 µl	50 25 5	1.00 % 1.50 % 3.00 %	0.50 % 1.00 % 2.00 %	721014	-
720340	12-channel	50-300 µl	5 µl	300 150 50	0.70 % 1.00 % 1.50 %	0.25 % 0.50 % 0.80 %	721014	-



Proline Single-Channel FIXED Volume, Performance Specifications and Ordering Details

Cat. No.	Channels	Volume Range (µl)	Increment (µl)	Test Volume (µl)	Inacc. (%)	Impr. (%)	Safe-Cone filters 50 pcs	
							Standard	Plus
722001	1-channel	5 µl	-	5	1.30 %	1.20 %	-	-
722004	1-channel	10 µl	-	10	0.80 %	0.80 %	-	-
722010	1-channel	20 µl	-	20	0.60 %	0.50 %	721008	721018
722015	1-channel	25 µl	-	25	0.50 %	0.30 %	721008	721018
722020	1-channel	50 µl	-	50	0.50 %	0.30 %	721008	721018
722025	1-channel	100 µl	-	100	0.50 %	0.30 %	721007	721017
722030	1-channel	200 µl	-	200	0.40 %	0.20 %	721007	721017
722035	1-channel	250 µl	-	250	0.40 %	0.20 %	721006	721016
722040	1-channel	500 µl	-	500	0.30 %	0.20 %	721006	721016
722045	1-channel	1000 µl	-	1000	0.30 %	0.20 %	721006	721016
722050	1-channel	2000 µl	-	2000	0.30 %	0.15 %	721005	721015
722055	1-channel	5000 µl	-	5000	0.30 %	0.15 %	721005	721015
721002	Biohit Carousel Stand for 5 pipettes							
725620	Biohit Linear Stand							

Due to manufacturing constraints we reserve the right to alter specifications without notice.

