



## Cardio-6T

Six Channel ECG

5.7" Color Touch Screen  
PC Review & Print  
JPG/PDF Reports  
Mass Storage Device  
Automatic Reporting System  
USB Flash Disk & SD Memory



The picture is only for reference, please make the object as the standard.

Principle : Korea

Assembly : China

### Product Features



Touch Mode



Thermal Paper



Quick Input



Unlimited Storage



USB & SD Drive



Color Display

### Specifications

1. 5.7" LCD screen.
2. ECG record modes.
3. Six Channel display modes.
4. 7 level printing thick can be chosen.
5. Calculate the heart rate automatically.
6. Thermal printer with 112 mm print paper.
7. With RS 232/USB communication interface.
8. Historical data can be reviewed and printed.
9. 12 lead ECG synchronous sampling system.
10. Adapt to 100-240V, 50/60Hz AC power supply.
11. With screen saver function when no work for a long time.
12. Auto report printing with comments (Normal/ Abnormal ECG).
13. Built-in SD memory card. Over 2000 patient's information can be stored.
14. Patient's information (ID number, Name, Gender, Age, Height, Weight, Blood pressure, etc.) can be input.
15. Automatic diagnosis function. ECG waveform can be measured and analyzed automatically.
16. Can upload real-time ECG data and history ECG data by connecting with ECG workstation. 12 lead synchronous printing can be realized.
17. Built-in rechargeable Li-ion battery 11.1V (3800mAh) with automatic protection function. ECG machine can work continuously for 10 hours and print over 1000 reports.

### Zhuhai Medical Science & Technology Co., Ltd.

Address: Room 910, B Area, No.8 Gangle Road, High-tech Zone, Zhuhai, China.  
Tel: 86-0756-6917978 Fax: 0756-6917976, Mobile: 18507562280