


Pelvic Floor Evaluation and Retraining



The Diversatek Healthcare Pelvic Floor Evaluation and Retraining System.

Diversatek Healthcare's comprehensive product portfolio for Pelvic Floor Evaluation and Retraining, including Anorectal Manometry, EMG Evaluation, and Pudendal Nerve Latency Testing. While there are many considerations involved in choosing the right products, know that the team at Diversatek Healthcare is ready to assist you every step of the way.


From Diversatek Healthcare's adaptive InSIGHT Ultima® System to BioVIEW® Acquisition and Analysis Software, including SyncVIEW™ that provides a multi-dimensional view of every maneuver, finding the right solution for your diagnostic and retraining needs has never been easier.




Anorectal Manometry

Diversatek Healthcare's anorectal manometry offering, including both directional and circumferential High Resolution models as well as the Clinical Innovations Air-Charged Catheters, allows the clinician to perform a broad spectrum of diagnostic tests, including:


- Push and Squeeze Maneuvers
- Recto-Anal Inhibitory Reflex (RAIR)
- Compliance Testing
- Sensation




Normal Push




Normal Squeeze




Squeeze Duration



Normal RAIR





Pelvic Floor Retraining

Functional disorders of the pelvic floor such as fecal incontinence and chronic constipation are disabling conditions which are increasing in prevalence. Most often the cause of these problems can be identified with a thorough physiological evaluation which includes Anorectal manometry, rectal sensation and compliance, pelvic floor EMG studies and Intra-rectal defecation studies. Once properly identified, these abnormalities can be treated.



One very useful tool is retraining the patient using various approaches based on desired outcome. The most important issue is having both the proper equipment and software to perform the task correctly. Diversatek Healthcare has developed solutions to enable accurate, user-friendly procedures for both patient and clinician.

Pudendal Nerve Latency Testing

Diversatek Healthcare offers the industry standard Sierra® Summit™ by Cadwell Laboratories, a complete system that measures the time relationship, or latency, between a stimulus of the Pudendal nerve and the external anal sphincter response.



Simplify Visualization

Diversatek Healthcare offers two styles of anorectal manometry high resolution probes. Choose a 23-channel, multi-directional probes or a 6 channel circumferential probe allowing visualization of the entire length of the anal canal.

BioVIEW® Software—with guided acquisition protocols and color pressure mapping—has simplified the complexity of information analysis and interpretation.

Innovations in Clinical Education

Diversatek University Online

Our online training platform contains free content on esophageal and anorectal manometric studies, as well as impedance/pH reflux monitoring studies. Included are tutorials providing step-by-step guidance to develop skills in data acquisition, study review and report generation. Simply go to **DiversatekHealthcare.com** to request log-in information.

Denver Training Center

Our Technical Research & Training Center specializes in clinical education and personalized training to make learning easy and convenient. Our comprehensive clinical instruction provides users with the knowledge and skills necessary to effectively acquire and analyze high resolution impedance manometry studies, impedance/pH reflux monitoring studies and anorectal manometry studies. Email us at **clinicaleducation@diversatekhc.com** or visit us online to learn more about our Denver course offerings.

Webinars

Diversatek Healthcare is proud to present a series of live, interactive discussions on topics related to esophageal function testing, impedance/pH reflux monitoring studies and anorectal manometry. Each webinar includes a didactic session followed by an open discussion. All webinars are recorded and posted to the Diversatek U online portal for easy reference. Access **DiversatekHealthcare.com** for upcoming webinar announcements.

The Diversatek Healthcare Review

Diversatek Healthcare Review, our monthly e-newsletter, features what's new at Diversatek University along with up-to-date product information. Every issue also includes our Clinical Insights, providing educational tips for Z/pH and HRiM analysis as well as answers to the most frequently asked questions.



Personalized Clinical Support

Onsite Training

Diversatek Healthcare Clinical Specialists deliver product support to suit your specific needs—on your schedule. Specialists are onsite at your facility to train and support you on your Diversatek Healthcare manometry or reflux monitoring equipment as you work through patient cases, acquire and analyze patient data, and create patient reports.

Virtual Coaching

Online and in real-time, Diversatek Healthcare Clinical Specialists work with you via screen sharing to provide study-specific data review and report generation coaching for your more difficult studies. Email us at **clinicalsupport@diversatekhc.com** to schedule a one-on-one session.

Diversatek[™] Healthcare

Advancing GI care by driving science, developing and delivering solutions, and providing unmatched clinical support.

DiversatekHealthcare.com

P 800.558.6408 or 414.265.7620

Corporate Headquarters

102 East Keefe Avenue
Milwaukee, WI 53212 USA

orders@diversatekhc.com
sales@diversatekhc.com

Technical Research & Training Center

9150 Commerce Center Circle, Suite 500
Highlands Ranch, CO 80129 USA

technicalsupport@diversatekhc.com
clinicalsupport@diversatekhc.com
clinicaleducation@diversatekhc.com

Learn more from
your personal product
representative at
800-558-6408