



DOMINUS

PRODUCT CATALOGUE



16-128 Slice CT System
Advanced Visualization Solution

A Better Insight into:

| Product quality

| Customer satisfaction

| Technical innovation




SG HealthCare is the global leader in the medical diagnostic equipment industry with world class technology advancements.

Today, we continue to spread the importance of "Seeing more" and "Knowing Better" as opposed to mediocre diagnosis. The information from our systems allows acute analysis, thus increasing the chance of saving lives.

X-ray, MRI & Ultrasound systems made by SG HealthCare are widely used all over the world from tiniest mobile clinics to national hospitals. SG HealthCare has proven itself to be revolutionizing the field of medical diagnostic systems and continues to help people all around the world see better with us.

Thank you.



DOMINUS

16/32/64S/64

Multi-function CT series















Bring more and more benefits to customers.



- High definition imaging(20 Lp/cm)with minimum dose.
- Easy to use with large 76 cm scan bore and 250kg patient couch.

▮ 32-row submillimeter detector with 20mm coverage realize fast scan and high quality images

▮ 5.3MHU/8.0MHU tube with 50kW/80kW generator provide long life and high profit for hospitals.



More Clarity with Less Dose



- 76cm large scan space provides greater openness for easier patient position.
- 64Rows/128 slice imaging platform generates full applications
- Fast workflow and excellent visualization solution
- Significant dose reduction and excellent clinical performance
- Cardiac application package

Optimize Your Financial Performance

Enhance Your Clinical Portfolio

Patient throughput optimization

Lower dose by design

DOMINUS

Built for clinical excellence by 32 slice platform







Enter the high-resolution CT world with the DOMINUS 32-slice capability and the clinical performance to enhance diagnostic confidence.

- 3.5M metal high capacity X-ray tube
- 350mA tube current selection
- 32 slice fast and precision scan
- Efficiency from positioning to post-processing.
- Perfect size of 0.6mm slice thickness detector
- Advance iterative reconstruction based on the noise model.
- 70kV advanced low dose scan

Advanced 128-slice CT designed to drive efficiency and full range application.

Top-level product with high configuration and economic profit to customer.






Vessel Analysis

One click Vessel Segmentation
And Measurement




Lung Nodule Analysis

Visualization Lung parenchyma
And segment nodule




Lung function analysis

Calculate the volume of emphysema,
Left lung, right lung and trachea




Tumor Assessment

Displaying tumor segmentation and
Measurement result



Body Perfusion

Providing whole-organ or single-location
Liver perfusion or tumor perfusion




Virtual Endoscopy


Fly-through for colon,
Trachea, vessel

Advanced visualization solution InsightVision post-processing workstation

The IVW is an advanced visualization.


Workstation that offers a integrated solution for multi-modalities and clinical domains to help you work quickly with increased diagnostic confidence





Brain Perfusion

Displaying the TDC, CBF, CBV, MTT, TTP images




CT-DSA

Display subtract result and Generate new data series




Dental Analysis

Creating panoramic image and sectional images




Cardiac Calcium Scoring

Measuring Calcium Score and Displaying Pseudo Color



Coronary Analysis

Vessel stenosis measurement; Coronary extraction and plaque analysis



Virtual Colonoscopy

Auto-segmentation Colon; Extraction Colon centerline














3D Volume Rendering

Image Cutting, Tissue Management, Volume Calculation




Full-rangeClinical
Applications






Advanced Visualization Solution

Making it easy for healthcare providers to share patient images, collaborate on clinical workflows, and communicate diagnosis report.



Post-processing collaboration, powered by Insight Vision clinical packages and tools.

Puts PACS system and workstation together with innovative cloud solution.



Global Service and Logistics Support

SG HealthCare delivers comprehensive service experience to our customers, while offering cost-effective and superior products.

- | Help our customers to control maintenance costs, mitigate medical risk, expand medical business and create more value for clinical diagnosis.
- | Remote diagnosis and predictive maintenance minimize downtime and maximize equipment life.



Corporate HQ

Add: 113 Jungdae-Ro, Songpa-Gu, Seoul, Korea

+82-2-449-6908

(Overseas Sales Dept. +82-70-4264-9704)

Sales@sghealthcare.com

service@sghealthcare.com